

Product Datasheet

ZNF423 antibody (orb387114)

Description	Rabbit polyclonal antibody to ZNF423
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ZNF423
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	145
Uniprot ID	Q2M1K9
Expiration Date	12 months from date of receipt.