

Product Datasheet

GEMIN6 antibody (orb387119)

Description	Rabbit polyclonal antibody to GEMIN6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human GEMIN6
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	19
Uniprot ID	Q8WXD5
Expiration Date	12 months from date of receipt.