

Product Datasheet

ZFP36 antibody (orb387123)

Description	Rabbit polyclonal antibody to ZFP36
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ZFP36
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	34
Uniprot ID	P26651
Expiration Date	12 months from date of receipt.