

Product Datasheet

IFT81 antibody (orb387128)

Description	Rabbit polyclonal antibody to IFT81
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human IFT81
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	80
Uniprot ID	Q8WYA0
Expiration Date	12 months from date of receipt.