

Product Datasheet

TNFRSF10D antibody (orb387148)

Description	Rabbit polyclonal antibody to TNFRSF10D
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human TNFRSF10D
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	42
Uniprot ID	Q9UBN6
Expiration Date	12 months from date of receipt.