

Product Datasheet

CD40LG antibody (orb387167)

Description	Rabbit polyclonal antibody to CD40LG
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CD40LG
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	29
Uniprot ID	P29965
Expiration Date	12 months from date of receipt.