

Product Datasheet

TRNP1 antibody (orb387182)

Description	Rabbit polyclonal antibody to TRNP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human TRNP1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	23
Uniprot ID	Q6NT89
Expiration Date	12 months from date of receipt.