

Product Datasheet

BAG4 antibody (orb387198)

Description	Rabbit polyclonal antibody to BAG4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human BAG4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	50
Uniprot ID	095429
Expiration Date	12 months from date of receipt.