

Product Datasheet

STAT5A antibody (orb387200)

Description	Rabbit polyclonal antibody to STAT5A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human STAT5A.
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	91
Uniprot ID	P42229
Expiration Date	12 months from date of receipt.