

Product Datasheet

KLF2 antibody (orb387206)

Description	Rabbit polyclonal antibody to KLF2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human KLF2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	37
Uniprot ID	Q9Y5W3
Expiration Date	12 months from date of receipt.