

Product Datasheet

SOX10 antibody (orb387213)

Description	Rabbit polyclonal antibody to SOX10
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SOX10
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	50
Uniprot ID	P56693
Expiration Date	12 months from date of receipt.