

Product Datasheet

ROBO1 antibody (orb387220)

Description	Rabbit polyclonal antibody to ROBO1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ROBO1
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	181
Uniprot ID	Q9Y6N7
Expiration Date	12 months from date of receipt.