

Product Datasheet

FZD2 antibody (orb387224)

Description	Rabbit polyclonal antibody to FZD2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human FZD2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	64
Uniprot ID	Q14332
Expiration Date	12 months from date of receipt.