

Product Datasheet

GITR antibody (orb387230)

Description	Rabbit polyclonal antibody to GITR
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse GITR
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	25
Uniprot ID	O35714
Expiration Date	12 months from date of receipt.