

Product Datasheet

ERCC1 antibody (orb387236)

Description	Rabbit polyclonal antibody to ERCC1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ERCC1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	33
Uniprot ID	P07992
Expiration Date	12 months from date of receipt.