

## Product Datasheet

### CHMP4C antibody (orb387259)

|                            |                                      |
|----------------------------|--------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CHMP4C |
| <b>Species/Host</b>        | Rabbit                               |
| <b>Reactivity</b>          | Human                                |
| <b>Conjugation</b>         | Unconjugated                         |
| <b>Tested Applications</b> | WB                                   |
| <b>Immunogen</b>           | N-terminal region of human CHMP4C    |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>               |
| <b>Storage</b>             | -20°C                                |
| <b>Note</b>                | For research use only                |
| <b>Clonality</b>           | Polyclonal                           |
| <b>Purity</b>              | Affinity purification                |
| <b>MW</b>                  | 26                                   |
| <b>Uniprot ID</b>          | <a href="#">Q96CF2</a>               |
| <b>Expiration Date</b>     | 12 months from date of receipt.      |