

Product Datasheet

MARF1 antibody (orb387260)

Description	Rabbit polyclonal antibody to MARF1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MARF1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	193
Uniprot ID	Q9Y4F3
Expiration Date	12 months from date of receipt.