

Product Datasheet

STX17 antibody (orb387273)

Description	Rabbit polyclonal antibody to STX17
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of mouse STX17
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	33
Uniprot ID	Q9D014
Expiration Date	12 months from date of receipt.