

Product Datasheet

SLX4 antibody (orb387278)

Description	Rabbit polyclonal antibody to SLX4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SLX4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	200
Uniprot ID	Q8IY92
Expiration Date	12 months from date of receipt.