

Product Datasheet

MFAP2 antibody (orb387284)

Description	Rabbit polyclonal antibody to MFAP2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human MFAP2
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	21
Uniprot ID	P55001
Expiration Date	12 months from date of receipt.