

Product Datasheet

MAP2K4 antibody (orb387292)

Description	Rabbit polyclonal antibody to MAP2K4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MAP2K4
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	44
Uniprot ID	P45985
Expiration Date	12 months from date of receipt.