

Product Datasheet

NR2C2 antibody (orb387303)

Description	Rabbit polyclonal antibody to NR2C2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human NR2C2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	65
Uniprot ID	P49116
Expiration Date	12 months from date of receipt.