


Product Datasheet

CST7 antibody (orb387304)

Description	Rabbit polyclonal antibody to CST7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CST7
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	16
Uniprot ID	076096
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot of mouse brain(lane 1), hum...