

## Product Datasheet

### FGF23 antibody (orb387306)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to FGF23 |
| <b>Species/Host</b>        | Mouse                              |
| <b>Reactivity</b>          | Human                              |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | ELISA, WB                          |
| <b>Immunogen</b>           | C-terminal region of human FGF23   |
| <b>Preservatives</b>       | 0.02% Thimerosal                   |
| <b>Storage</b>             | -20 °C                             |
| <b>Note</b>                | For research use only              |
| <b>Clonality</b>           | Monoclonal                         |
| <b>Purity</b>              | Affinity purification              |
| <b>MW</b>                  | 28                                 |
| <b>Uniprot ID</b>          | <b>Q9GZV9</b>                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |