

Product Datasheet

MEAF6 antibody (orb387329)

Description	Rabbit polyclonal antibody to MEAF6
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MEAF6
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	22
Uniprot ID	Q9HAF1
Expiration Date	12 months from date of receipt.