

Product Datasheet

CALB1 antibody (orb387338)

Description	Rabbit polyclonal antibody to CALB1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CALB1
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	30
Uniprot ID	P05937
Expiration Date	12 months from date of receipt.