

Product Datasheet

P4HB antibody (orb387349)

Description	Rabbit polyclonal antibody to P4HB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human P4HB
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	57
Uniprot ID	P07237
Expiration Date	12 months from date of receipt.