

Product Datasheet

PCSK9 antibody (orb387350)

Description	Rabbit polyclonal antibody to PCSK9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PCSK9
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	74
Uniprot ID	Q8NBP7
Expiration Date	12 months from date of receipt.