

Product Datasheet

HSD17B2 antibody (orb387351)

Description	Rabbit polyclonal antibody to HSD17B2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	various regions of human HSD17B2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	43
Uniprot ID	P37059
Expiration Date	12 months from date of receipt.