

Product Datasheet

PCGF1 antibody (orb387355)

Description	Rabbit polyclonal antibody to PCGF1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human PCGF1
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	30
Uniprot ID	Q9BSM1
Expiration Date	12 months from date of receipt.