

Product Datasheet

PAIP2 antibody (orb387357)

Description	Rabbit polyclonal antibody to PAIP2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human PAIP2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	15
Uniprot ID	Q9BPZ3
Expiration Date	12 months from date of receipt.