

Product Datasheet

TPSAB1 antibody (orb387358)

Description	Rabbit polyclonal antibody to TPSAB1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TPSAB1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	30
Uniprot ID	Q15661
Expiration Date	12 months from date of receipt.