


Product Datasheet

BEX1 antibody (orb387375)

Description	Brain-expressed X-linked protein 1 (BEX1) is a signaling adapter molecule involved in p75NTR/NGFR si...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human BEX1
Preservatives	0.01% NaN3
Form/Appearance	Supplied in 1 x PBS (pH 7.4), 100 ug/ml rAlbumin, 40% Glycerol, 0.01% NaN3
Storage	Store at -20 °C. Stable for 6 months from the date of receipt
Note	For research use only
Clonality	Polyclonal
Purity	Antibodies were purified by affinity purification using immunogen
MW	15 kDa
Uniprot ID	Q9HBH7
Dilution Range	WB: 1:1,000-1:3,000 (detect endogenous protein*)
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...