

Product Datasheet

NHP2 antibody (orb387377)

Description	Rabbit polyclonal antibody to NHP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NHP2
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	17
Uniprot ID	Q9NX24
Expiration Date	12 months from date of receipt.