

Product Datasheet

SPTAN1 antibody (orb387383)

Description	Rabbit polyclonal antibody to SPTAN1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SPTAN1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	285
Uniprot ID	Q13813
Expiration Date	12 months from date of receipt.