

Product Datasheet

ZNRF2 antibody (orb387388)

Description	Rabbit polyclonal antibody to ZNRF2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ZNRF2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	24
Uniprot ID	Q8NHG8
Expiration Date	12 months from date of receipt.