

Product Datasheet

Tubulin A antibody (orb387390)

Description	Mouse monoclonal antibody to Tubulin A
Species/Host	Mouse
Reactivity	Human, Mouse, Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human Tubulin A1A
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Monoclonal
Purity	Affinity purification
MW	50
Uniprot ID	Q71U36
Expiration Date	12 months from date of receipt.