

Product Datasheet

CASC5 antibody (orb387395)

Description	Rabbit polyclonal antibody to CASC5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CASC5
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	265
Uniprot ID	Q8NG31
Expiration Date	12 months from date of receipt.