


Product Datasheet

PAK4 antibody (orb38740)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	PAK4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64072
Uniprot ID	O96013
NCBI	NP_001014832.1 , NP_001014834.1 , NP_001014835.1 , NP_005875.1 , NP_001014831.1
Dilution Range	IHC-P: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of in mouse small intestine tissue.

