

Product Datasheet

BMP11 antibody (orb387420)

Description	Rabbit polyclonal antibody to BMP11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse BMP11
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	45
Uniprot ID	Q9Z1W4
Expiration Date	12 months from date of receipt.