

Product Datasheet

S1PR4 antibody (orb387438)

Description	Rabbit polyclonal antibody to S1PR4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human S1PR4
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	42
Uniprot ID	095977
Expiration Date	12 months from date of receipt.