

Product Datasheet

RPS6KA4 antibody (orb387439)

Description	Rabbit polyclonal antibody to RPS6KA4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human RPS6KA4
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	86
Uniprot ID	075676
Expiration Date	12 months from date of receipt.