



Product Datasheet

PIM1 antibody (orb38744)

Description	Rabbit polyclonal antibody to PIM1 or pim-1 whose expression can be detected in the blood cells but ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	PIM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35686
Uniprot ID	P11309
NCBI	NP_002639.1 , NP_001230115.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...