

Product Datasheet

SMAD1 antibody (orb387440)

Description	Rabbit polyclonal antibody to SMAD1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SMAD1
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	52
Uniprot ID	Q14797
Expiration Date	12 months from date of receipt.