




Product Datasheet

PIM2 antibody (orb38745)

Description	Rabbit polyclonal antibody to PIM2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Immunogen	Synthetic Peptide
Target	PIM2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34190
Uniprot ID	Q9P1W9
NCBI	NP_006866.2
Dilution Range	WB: 1:1000, IP: 1:50-100, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of MV4;11 C...



Western blot analysis of in HeLa cell li...