

Product Datasheet

PAIP1 antibody (orb387464)

Description	Rabbit polyclonal antibody to PAIP1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PAIP1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	54
Uniprot ID	Q9H074
Expiration Date	12 months from date of receipt.