

Product Datasheet

MYO5C antibody (orb387469)

Description	Rabbit polyclonal antibody to MYO5C
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MYO5C
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	203
Uniprot ID	Q9NQX4
Expiration Date	12 months from date of receipt.