

Product Datasheet

TRPS1 antibody (orb387471)

Description	Rabbit polyclonal antibody to TRPS1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human TRPS1
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	142
Uniprot ID	Q9UHF7
Expiration Date	12 months from date of receipt.