

Product Datasheet

SETD1B antibody (orb387480)

Description	Rabbit polyclonal antibody to SETD1B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse SETD1B.
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	215
Uniprot ID	Q8CFT2
Expiration Date	12 months from date of receipt.