

Product Datasheet

L3MBTL antibody (orb387502)

Description	Rabbit polyclonal antibody to L3MBTL
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human L3MBTL
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	84
Uniprot ID	Q9Y468
Expiration Date	12 months from date of receipt.